



[Articles](#) » [Measurement of Viscosity in a Vertical Falling Ball Viscometer](#)

## Measurement of Viscosity in a Vertical Falling Ball Viscometer

Ping Yuan

Ben-Yuan Lin


Posted: October 27, 2008

### Experimental

The experimental devices were arranged as shown in *Figure 2a*, which includes the fall tube containing the liquid to be measured, a tubular jacket for thermal control, a water reservoir with a heater, a temperature control system, a pump for circulation, and a camera. The camera has four million pixels, and captures an image per 1/30 sec. The temperature control system can control the water temperature deviation to less than 0.5 °C via feedback from a sensor in the tubular jacket.

Figure 2b is a photograph of the entire device except for the camera. In this figure, the device on the right is the water reservoir with a heater and a pump, which connects to the tubular jacket via soft pipes. The device in the center of Figure 2b is the power and control box, which supplies the power to the pump and heater in the reservoir and controls the on/off of the heater to maintain a set temperature of water. The device on the left is the main structure for the falling ball. The fall tube and the tubular jacket can be inverted by rotating them together about their mounting point on the stand in order to return the ball to the starting position. On the top of the fall tube is a pin to prevent the ball from falling, and it can be released by the control on the control box. The stand is equipped with a level and feet with leveling screws for adjusting the fall tube so that it is vertical to the ground.

The experimental procedure is as follows: 1) Pour the liquid under test into the fall tube and place the ball into the tube. Then, tighten the screw on the top of the fall tube. 2) Set the measurement temperature and turn on the heater and pump. It takes at least 15 min for the tubular jacket temperature to reach the specified test temperature.



### Claim your Complimentary Subscription

If you are a scientist, you are eligible for a full subscription.

[Click here to get started »](#)



### Related Products



#### DV-I Prime Continuous Sensing Digital Viscometer

The DV-I Prime Viscometer has continuous sensing capability for rapid viscosity measurement makes this economical digital viscometer a standout. ...



#### Falling Ball Viscometer

The Höppler principle is used to measure the viscosity of Newtonian liquid by measuring the time required for a ball to fall under gravity through a ...

### Related Product Categories

[Falling Ball Viscometer / Falling Sphere Viscometers »](#)

[Rotational Viscometer / Rotary Viscometers »](#)

### Subscribe to eNewsletters and Email Alerts

Keep up with our latest articles, news and events. Plus, get special offers and more delivered to your inbox.

Your Email Address

3) Rotate the fall tube 180° and fix the pin for holding the falling ball. Then, rotate it back and ready to fall ball. 4) Set up the camera to be perpendicular to the fall tube, and begin to record. 5) Release the pin for dropping the ball, and keep the camera record until the falling ball passes through the marked end line of the fall tube. 6) Record the falling time from the beginning to the end, and follow the calculation step in the above section to obtain the dynamic viscosity.



Results and discussion

Table 1 - Measuring and calculating results of SAE 30 at 40 °C in the experiment

Exp. No.	Falling length (5cm)			Falling length (10cm)		
	time(s)	Re	$\mu$	time(s)	Re	$\mu$
1	1.23	3.3991	0.1044	2.40	3.5170	0.1034
2	1.17	3.7708	0.0989	2.27	3.9395	0.0976
3	1.20	3.5776	0.1016	2.33	3.7355	0.1003
4	1.17	3.7708	0.0989	2.27	3.9395	0.0976
5	1.13	4.0540	0.0952	2.23	4.0850	0.0958
6	1.23	3.3991	0.1044	2.40	3.5170	0.1034
7	1.17	3.7708	0.0989	2.30	3.8355	0.0989
8	1.20	3.5776	0.1016	2.37	3.6082	0.1020
9	1.20	3.5776	0.1016	2.37	3.6082	0.1020
10	1.20	3.5776	0.1016	2.37	3.6082	0.1020
Exp. No.	Falling length (15cm)			Falling length (20cm)		
	time(s)	Re	$\mu$	time(s)	Re	$\mu$
1	3.57	3.7681	0.1001	4.67	3.6869	0.1014
2	3.40	3.9268	0.0980	4.47	4.0276	0.0969
3	3.43	3.8575	0.0989	4.47	4.0276	0.0969
4	3.37	3.9979	0.0971	4.43	4.1014	0.0961
5	3.37	3.9979	0.0971	4.40	4.1581	0.0954
6	3.60	3.4977	0.1039	4.57	3.8516	0.0991
7	3.47	3.7680	0.1001	4.50	3.9735	0.0976
8	3.57	3.7681	0.1001	4.60	3.8010	0.0998
9	3.57	3.7681	0.1001	4.67	3.6869	0.1014
10	3.63	3.4395	0.1048	4.70	3.6396	0.1020

The present study uses an oil of SAE 30 as the sample liquid for measurement. It belongs to the CD/SF level in the classification of the American Petroleum Institute (API), and its kinetic viscosity is 111.2 cSt at 40 °C, which is equal to 0.0973 Pa·s of dynamic viscosity. The sample oil has a density of 894.4 kg/m<sup>3</sup> through measurement by an electrical balance with a precision of 1 mg. Moreover, the study utilizes a plastic sphere ball as the falling ball, which has a diameter of 9.756 mm, out-of-roundness of 0.09 mm, and density of 1010 kg/m<sup>3</sup>. The experiment proceeds in a 40 °C tubular jacket temperature, with four different falling distances of 5, 10, 15, and 20 cm. Furthermore, this study was repeated 10 times for measuring the fall time in each case; the experimental results are listed in Table 1.




Figure 3 - Relative error of measuring dynamic viscosity to the reference value in different falling distance.

this experiment.<sup>10</sup>Table 2 lists the mean of results; experimental standard deviation; uncertainty; and confidence internal, which are defined as follows:

Table 2 - Uncertainty and the confidence internal in this experiment

Fall Distance	$\bar{\mu}$	$s$	$u$	$\pm 4.09u$
5cm	0.10071	0.28006E-4	0.88562E-3	$\pm 0.00362$
10cm	0.10030	0.26843E-4	0.84885E-3	$\pm 0.00347$
15cm	0.10002	0.25890E-4	0.81900E-3	$\pm 0.00335$
20cm	0.09866	0.24114E-4	0.76255E-3	$\pm 0.00312$

$$\bar{\mu} = \frac{1}{n} \sum_{j=1}^n \mu_j \tag{10}$$

Find us Follow us

advertisement

**NEW!**

**labcompare**

**Clinical Diagnostics Assay Kits**

Hundreds of available assays from several suppliers.

**Request Quote and Compare Side-by-Side:**

- Sample
- Assay Method
- Detection
- Sensitivity
- And more!

**Make your research less complicated...Search Now!**

Product Guide

- Analytical Instruments
- Chemical Analysis Equipment
- Clinical Diagnostics
- Environmental Test Equipment
- Food Testing Equipment
- Forensic Laboratory Equipment
- General Laboratory Equipment
- Lasers and Optical Components
- Mass Spectrometry
- Materials Research
- Microscopy
- Pharmaceutical Lab Equipment
- Spectroscopy
- Test and Measurement

more products from labcompare »

Industry News

Saint-Gobain Seals Group Exhibiting at Pittcon 2014 to Showcase Their Innovative Sealing and Polymer Solutions in Life Sciences Industry

Wednesday, February 19, 2014

Eppendorf Partners with Biomanufacturing Education and Training Center at WPI

Wednesday, February 19, 2014

AFM-IR Used by French Researchers to Identify Best Microbes for Biofuel Production

Wednesday, February 19, 2014

MATHESON Acquires Continental Carbonic Products, Inc.

Wednesday, February 19, 2014

more news »

Most Viewed Articles

Most Recent Articles

Determination of Palladium in Water and Geological Samples by Flame Atomic Absorption Spectrophotometry

Thursday, November 28, 2013

$$s = \sqrt{\frac{1}{n-1} \sum_{j=1}^n (\mu_j - \bar{\mu})^2} \quad (11)$$

$$u = \sqrt{\frac{1}{n(n-1)} \sum_{j=1}^n (\mu_j - \bar{\mu})^2} \quad (12)$$

$$\mu = \bar{\mu} \pm 4.09u \quad (13)$$

In Table 2, the uncertainty and confidence interval decrease with the increase in falling distance, and the confidence interval is  $\pm 0.0034$  with a confidence level of 99.73%. This demonstrates that the error of measurement is about 6%, which is an acceptable experiment error.


## Comments

4 Comments American Laboratory  Login ▾


Sort by Best ▾ Share  Favorite ★




Join the discussion...




Sina Changizian • 2 years ago  
good  
tnx alot!  
1 ^ | ▾ • Reply • Share ›



Zeeshanhaider14 • 2 years ago  
good  
1 ^ | ▾ • Reply • Share ›



Omar Rifaie • 4 months ago  
Perfect, sehr schön. :)  
Thx a lot  
^ | ▾ • Reply • Share ›




Marwan Taher ➔ Omar Rifaie • 4 months ago  
Fluid mechanics hausaufgaben? ;)  
^ | ▾ • Reply • Share ›

### ALSO ON AMERICAN LABORATORY

#### Three Reasons Not to Throw Out Your LIMS

1 comment • 10 months ago



LIMSExpert.com — So here's the thing, the statement -- "this custom-coding (programming) is the bane of the LIMS ...

#### April 2013

14 comments • 11 months ago




Irene — Escherichia coli

### WHAT'S THIS?

#### Spectroscopic Ellipsometry Characterization of Thin Films Used ...


1 comment • 8 months ago



food packaging — First of all i would like to thank you for the great and informative entry.

#### Field-Flow Fractionation Supporting Consumer Safety Evaluation of Silver ...

1 comment • 11 months ago



Verbsat — A complete bibliography of every (known) paper published on FFF may be found at Wyatt.com.

### New Technologies in Human Genetic Analysis

Wednesday, April 04, 2007

### Measurement of Viscosity in a Vertical Falling Ball Viscometer

Monday, October 27, 2008

### Positional Fluctuation of IR Absorption Peaks: Frequency Shift of a Single Band or Relative Intensity Changes of Overlapped Bands?

Thursday, March 10, 2011

### Determination of Association ( $k_{on}$ ) and Dissociation ( $k_{off}$ ) Rates of Spiperone on the Dopamine D2 Receptor Using a Platform for GPCR Applications

Saturday, October 01, 2011

## Write for Us

American Laboratory wants to hear from you. Submit your case studies, applications, laboratory management tips, editorial, or other manuscripts.

[Learn more about writing for us](#)